

Step 1: Upload Images

- 📁 **Supported Types:** B-mode images and CEUS video *PNG/JPG/MP4*
 - 🖼️ You can upload multiple images at once, drag & drop supported.
 - 🖋️ Please make sure image resolution is **400×600** for best results.
 - 🔍 Preview images before submission to avoid input mistakes.
 - 🆔 File names must be anonymized *nopatientinfo*.
 - ⌚ CEUS video should be at least **10 seconds** long to capture perfusion.
 - 📌 *Tip:* The system auto-selects key CEUS frames for optimal performance.
-

Step 2: Run AI Inference

- 🤖 Click **Run Diagnosis** to launch the model.
 - 🧠 The model performs tumor segmentation and classifies lesion as **benign or malignant**.
 - 🏷️ Extracted features include:
 - 📈 Time-intensity curve
 - 📊 Peak intensity & time-to-peak
 - 📍 Tumor center & size
 - ⚡ Inference takes ~5s per case and runs on either local or cloud backend.
 - 🔍 Results will be visualized immediately after processing completes.
-

Step 3: Generate Diagnostic Report

- 📄 All results are compiled into a structured, printable report.
- 🖼️ Includes original image with overlay, tumor location, and AI classification.
- 📁 One-click export available in PDF and JSON formats.
- 🗄️ Each report has a unique case ID and is automatically saved in the case archive.
- 🔒 Reports can be securely shared with medical teams or uploaded to patient records.
- 🛠️ *Tip:* You can customize the report format in the **Settings** panel.